

Improved thermophilic genome editing tools for engineering *Clostridium thermocellum*

Background

- Discovery and validation of functional thermophilic genetic tools are needed to accelerate genome editing in the cellulolytic chassis bacterium *Clostridium thermocellum*. Standard tools do not function at thermophilic temperatures.

Approach

- Characterized the native Type I-B CRISPR-Cas system as well as evaluated thermophilic Type II systems to enable CRISPR-based genome editing
- Tested potential thermophilic homologs of Lambda Red and RecE/T recombineering systems to improve homologous recombination efficiencies



Outcome

- Identified functional protospacer adjacent motif (PAM) sequences for efficient targeting for cutting by the Type I-B system
- A Type II system previously isolated from *Geobacillus* was found to be highly effective for genome targeting and cutting
- Lambda Red beta/exo homologs from *Acidithiobacillus caldus* improved homologous recombination efficiencies by 35-fold
- Developed a markerless, 2-step system for generating mutants with up to 94% efficiency

Significance

- These tools allow for accelerated specifically targeted mutant generation (knockouts, knockins, SNPs) in *C. thermocellum* and potentially other thermophiles

2-step CRISPR-based Genome Editing



*Correct transformants = TE * % edited colonies

** Cas9 nickase variant was used

*** Rec = *A. caldus* Exo/Beta recombineering machinery